




Presents:

HYPERBARIC OXYGEN THERAPY


Examples of HBOT Enhancing Treatment & Care

Severe Effusive Canine Pancreatitis
Associated with Progressive Pleural Fluid
Treated for 14 days with standard, conservative therapy
prior to referral for Hyperbaric Oxygen Therapy

Pre Hyperbaric Oxygen Therapy



Post Hyperbaric Oxygen Therapy



© RKLyman, LLC. All Rights Reserved.

In addition to standard supportive therapy, Hyperbaric Oxygen Therapy (HBOT) was given at 1.5 ATA for 1 hour on:

- Day 1 twice daily**
- Day 2 twice daily**
- Day 3 twice daily**
- Day 4 twice daily**

For more information about this HBOT case and many more indications, please [click here](#) and you will be promptly contacted by Dr. Ronald Lyman, DVM, Diplomate ACVIM, and Medical Director of Veterinary Hyperbaric Safety.